



1958 / 2008 YEARS 

## 3C86/3C96 Compact™ Microprocessor Knee

New Options – Same Remarkable Stability

There's more than meets the eye with the new Compact Microprocessor Knee. For patients who meet Mobility Grades 2 or 3, Compact offers features such as unprecedented stability and a remote control for Static Stance to lock the knee for extended standing.

*OAO Bock*®

QUALITY FOR LIFE

# 3C86/3C96

## Compact™ Microprocessor Knee



### Features and Benefits

- Static Stance to lock the knee between 0° and 30°
- Remote control to activate and deactivate Static Stance
- Vibration alerts the user when Static Stance is engaged/disengaged
- Option of using a clip to carry the remote control or using a low-profile mounting plate to attach it to the socket
- Microprocessor control for superb stability
- Easier initiation of swing phase

### Patient Indications

The Compact is designed for amputees who meet Mobility Grades 2 or 3 and weigh up to 275 lbs (125 kg). This knee is designed for people who are (or have the potential to be) community ambulators, need a high degree of stance phase stability, but would not benefit from highly variable swing phase dynamics. The Static Stance feature makes the Compact an excellent option for people who need to stand for extended periods of time.

### Ordering Information

For more information about the Compact or other Otto Bock HealthCare products or to place an order, please call your local Sales Representative at 800.328.4058 (U.S.) or 800.665.3327 (Canada). You can also find more information online at [www.ottobockus.com](http://www.ottobockus.com) (U.S.) or [www.ottobock.ca](http://www.ottobock.ca) (Canada).

### Warranty and Service

36-month warranty. Service 24 months after purchase. Please refer to the Instructions for Use for complete information.

