



QUALITY FOR LIFE

C-Soft 2.2 System Requirements

- C-Soft 2.2 requires Windows 7 (32 or 64 bit), Windows Vista (32 or 64 bit) or Windows XP (SP2 or SP3)
- Netbooks **must have** a minimum screen resolution of 1024 x 768. C-Soft will not display correctly at a lower resolution.

Minimum Computer Requirements

- A PC with a 1000 MHz Pentium III processor
- 512 MB RAM
- Screen resolution of 1024 x 768 or higher
- 1 GB of free hard drive space
- CD-ROM drive
- RS232 (Serial) or USB connection

C-Soft 2.0 System Requirements

- C-Soft 2.0 **is not** compatible with 64-bit operating systems. Upgrade to C-Soft 2.2.
- C-Soft 2.0 requires Windows Vista **32-bit** or Windows XP (SP2 or SP3).
- Netbooks **must have** a minimum screen resolution of 1024 x 768. C-Soft will not display correctly at a lower resolution.

Minimum Computer Requirements

See C-Soft 2.2 System Requirements

C-Soft 1.0 and C-Soft 1.2 System Requirements

- C-Soft 1.0 and C-Soft 1.2 are **NOT** compatible with 64-bit operating systems.
- C-Soft 1.0 requires Windows XP SP2 or SP3.
- C-Soft 1.2 requires Windows Vista 32-bit or Windows XP (SP2 or SP3). If you have purchased C-Soft 1.0 and need to upgrade to C-Soft 1.2, contact Vicki Wermager at 800-328-4058 x5172.
- Netbooks **must have** a minimum screen resolution of 1024x768. C-Soft will not display correctly at a lower resolution.

Minimum Computer Requirements

- A PC with a 400 MHz Pentium II processor
- 128 MB RAM
- Screen resolution of 1024 x 768 or higher
- 120 MB of free hard drive space
- CD-ROM drive
- RS232 (Serial) or USB connection

If you have any questions regarding C-Soft system requirements, please call Otto Bock @ 800-328-4058.