

# C-Soft System Requirements

## Supported Operating Systems and Minimum System Requirements

---

This document contains the minimum system requirements and supported operating systems for C-Soft 2.4, C-Soft 2.2, C-Soft 2.0, C-Soft 1.2, and C-Soft 1.0.

**Note:** *Otto Bock Software is not compatible with Macintosh computers and should not be run on a Macintosh computer using PC emulator software.*

## C-Soft 2.4/2.2 System Requirements

### Supported Operating Systems (32-bit)

- Windows XP Professional (SP2 or later)
- Windows XP Home Edition (SP2 or later)
- Vista, all editions except Vista Starter
- Windows 7, all editions except Windows 7 Starter

### Supported Operating Systems (64-bit)

- Windows 7, all editions except Windows 7 Starter

### Recommended Computer Requirements

For 32-bit platforms:

- PC with a Pentium IV 650 / 3.4 GHz processor (x86) or faster
- At least 1 GB RAM

For 64-bit platforms:

- PC with an Intel Core2 (x64) or faster
- At least 2 GB RAM

For 32-bit and 64-bit platforms:

- At least 1 GB free hard drive space
- Screen resolution of 1024 x 768 at 96 DPI
- CD-ROM drive
- Mouse and keyboard

### Minimum Computer Requirements

- PC with a Pentium III /1 GHz processor
- 512 MB RAM
- Screen resolution of 1024 x 768 at 96 DPI
- 1 GB free hard drive space
- CD-ROM drive
- Mouse and keyboard

## C-Soft 2.0 System Requirements

C-Soft 2.0 is not compatible with Windows 7. If you are running Windows 7 and you have C-Soft 2.0 you have two options to upgrade:

- Install C-Soft 2.0 on Windows 7 and then install the 2.2 update located on the Otto Bock Software Support website  
<http://academy.ottobockus.com/software/software.aspx>.
- Call customer service and order C-Soft 2.2 from Otto Bock at no charge

### **Supported Operating Systems (32-bit)**

- Windows XP Professional (SP2 or later)
- Windows XP Home Edition (SP2 or later)
- Vista, all editions except Vista Starter

### **Minimum Computer Requirements**

- PC with a Pentium III /1 GHz processor
- 512 MB RAM
- Screen resolution of 1024 x 768 at 96 DPI
- 1 GB free hard drive space
- CD-ROM drive
- Mouse and keyboard

## C-Soft 1.2 System Requirements

### **Supported Operating Systems (32-bit)**

- Windows XP Professional (SP2 or later)
- Windows XP Home Edition (SP2 or later)
- Vista, all editions except Vista Starter

### **Minimum Computer Requirements**

- A PC with a 400 MHz Pentium II processor
- 128 MB RAM
- Screen resolution of 1024 x 768 or higher
- 120 MB of free hard drive space
- CD-ROM drive
- RS232 (Serial) or USB connection

## C-Soft 1.0 System Requirements

### **Supported Operating Systems (32-bit)**

- Windows XP Professional (SP2 or later)
- Windows XP Home Edition (SP2 or later)

### **Minimum Computer Requirements**

- A PC with a 400 MHz Pentium II processor
- 128 MB RAM
- Screen resolution of 1024 x 768 or higher
- 120 MB of free hard drive space
- CD-ROM drive
- RS232 (Serial) or USB connection