

ElbowSoft System Requirements

Supported Operating Systems and Minimum System Requirements

This document contains the minimum system requirements and supported operating systems for ElbowSoft 1.6.

Note: *Otto Bock Software is not compatible with Macintosh computers and should not be run on a Macintosh computer using PC emulator software.*

ElbowSoft 1.6 System Requirements

Supported Operating Systems (32-bit)

- Windows XP Professional (SP2 or later)
- Windows XP Home Edition (SP2 or later)
- Vista, all editions except Vista Starter
- Windows 7, all editions except Windows 7 Starter

Supported Operating Systems (64-bit)

- Windows 7, all editions except Windows 7 Starter

Recommended Computer Requirements

For 32-bit platforms:

- PC with a Pentium IV 650 / 3.4 GHz processor (x86) or faster
- At least 1 GB RAM

For 64-bit platforms:

- PC with an Intel Core2 (x64) or faster
- At least 2 GB RAM

For 32-bit and 64-bit platforms:

- At least 1 GB free hard drive space
- Screen resolution of 1024 x 768 at 96 DPI
- CD-ROM drive

Minimum Computer Requirements

- PC with a Pentium III /1 GHz processor
- 512 MB RAM
- Screen resolution of 1024 x 768 at 96 DPI
- 1 GB free hard drive space
- CD-ROM drive