

PAULA System Requirements

Supported Operating Systems and Minimum System Requirements

This document contains the minimum system requirements and supported operating systems for PAULA 1.0.

Note: *Otto Bock Software is not compatible with Macintosh computers and should not be run on a Macintosh computer using PC emulator software.*

PAULA 1.0 System Requirements

Supported Operating Systems (32-bit)

- Windows XP Professional (SP2 or later)
- Windows XP Home Edition (SP2 or later)
- Vista, all editions except Vista Starter

Supported Operating Systems (64-bit)

- None

Recommended Computer Requirements

- PC with a Pentium IV 650 / 3.4 GHz or faster processor
- 1 GB RAM
- At least 1 GB free hard drive space
- Screen resolution of 1024 x 768 at 96 DPI
- CD-ROM drive
- Mouse and keyboard

Minimum Computer Requirements

- PC with a Pentium III /1 GHz processor
- 512 MB RAM
- 1 GB free hard drive space
- Screen resolution of 1024 x 768 at 96 DPI
- CD-ROM drive
- Mouse and keyboard

If you have any questions regarding C-Soft system requirements, please call Otto Bock @ 800-328-4058.