



**QUALITY FOR LIFE**

## OBSS 6.5 System Requirements ISOTRAK / Patriot System

**NOTE:** The ISOTRAK digitizing equipment has been discontinued and replaced with the Patriot system, the equipment is not interchangeable.

The OBSS can use either system although **the Patriot System will only work with Windows 2000, Windows XP or Windows Vista**. The ISOTRAK system will work on Windows 98SE, Windows ME, Windows 2000, Windows XP or Windows Vista. If you have any questions regarding the system requirements, please call Otto Bock at 800-328-4058.

### *Recommended Computer Requirements*

- A PC with 1.6 GHz (or faster) Pentium IV processor
- 256 MB RAM (or more)
- 120 MB or more of free hard drive space
- Graphics card (NVIDIA or ATI) with at least 16 MB video memory
- Connection to internet with email
- Free serial or USB port

### *Minimum Computer Requirements*

- A PC with 350 MHz Pentium II processor
- 128 MB RAM
- 120 MB or more of free hard drive space
- Graphics card with at least 16 MB video memory
- Connection to internet with email
- Free serial or USB port

