



**QUALITY FOR LIFE**

## TF / TT System Requirements

TF and TT design requires one of the following operating systems: Windows 2000, Windows XP or Vista. If you have any questions regarding the system requirements please call Otto Bock at 800-328-4058.

### *Recommended Computer Requirements*

- A PC with Pentium IV 650 / 3.4 GHz-Processor
- 1024 MB RAM (or more)
- 1 GB or more of free hard disk storage space
- CD-ROM or DVD-ROM drive
- Minimum resolution of 1024 x 768 pixels
- Email and internet provider

### *Minimum Computer Requirements*

- A PC with 1000 MHz Pentium III processor
- 512 MB RAM
- 1 GB or more of free hard disk storage space
- CD-ROM drive
- Graphics card with support for OpenGL 3D hardware acceleration  
Resolution of 1024 x 768 pixels or higher

