



QUALITY FOR LIFE

TT Design Supported Cameras

The Otto Bock TT Design software only supports cameras listed in this document. Although Canon manufactures numerous camera models only the cameras listed have been calibrated to work with TT design.

Supported Cameras

Camera	Date Added	TT Version
Kodak DC260, DC265, DC290	Before 12/8/2002	-
Canon A20	Before 12/8/2002	-
Canon A40, A60	3/9/2003	3.1.5
Canon A75	3/26/2004	3.9.3
Canon A85	5/20/2005	4.0.16
Canon IXUS60	8/10/2006	5
Canon PowerShot SD1100 IS (IXUS80)	8/19/2008	5.14
Canon PowerShot SD1200 IS (IXUS95)	8/1/2010	5.18

If you need additional technical support, contact Otto Bock at 800-328-4058.